



DIGITAL ULTRASONIC CLEANER

LDUC-A17

Digital Ultrasonic Cleaner LDUC-A17

Digital Ultrasonic Cleaner LDUC-A17 is a compact and robust 22 L tanked structure that creates waves in the tank water to clean glassware and instruments from grime. The inner steel chamber is resistant to cavitation bubbles and strong solvents when intense cleaning is necessary. Capable of adjusting the power exerted, the ultrasonic cleaner also protects fragile equipment. Equipped with degassing and degreasing function for removal of dissolved gasses and cleansing of the tank.

Features

- ❖ 22 L capacity tank
- ❖ Integrated digital circuit for timer and heater
- ❖ Temperature range : 0 °C ~ 80 °C
- ❖ Cleaning duration : 1 ~ 99 mins
- ❖ Gapless purge tank
- ❖ One key access for auto stand-by , sleep and wake up mode
- ❖ Efficient cleaning can be done by distilled water
- ❖ Precise digital processed temperature controls
- ❖ Optimum durability under prolonged and heavy usage

Application

Ideal for cleaning lenses and other optical parts, dental and surgical instruments, electronics, tools, pens, firearms, jewelry, and golf clubs. It is also used in manufacturing and maintenance of aerospace, engineering, filtration, orthopedic implants, pharmaceutical, plating and surface finishing, printing technology and equipment etc.

Specifications

| Model No. | LDUC-A17 |
|----------------------------|--------------------------|
| Capacity | 22 L |
| Tank size (L x W x H) | 500 mm x 300 mm x 150 mm |
| Machine size (L x W x H) | 500 mm x 330 mm x 330 mm |
| Ultrasonic power | 480 W |
| Heater power | 500 W |
| Temperature set | 0 °C ~ 80 °C |
| Digital timer | 1 - 99 mins |
| Control mode | Digital |
| Construction material | Stainless Steel |
| Voltage | AC 220 V / 50 Hz |

Optional Accessories

| | |
|-----------------|-----|
| Cleaning Basket | Yes |
| Drainage | Yes |