



Rapid Low Temperature Circulating Bath

LRTB-A40

Rapid Low Temperature Circulating Bath LRTB-A40

Rapid Low Temperature Circulating Bath LRTB-A40 is a microcomputer controlled 7 L bath with refrigeration system and a temperature range of -40 °C to 200 °C. Equipped with dual red and green LED windows with upper window measuring red value and lower measuring green value. Features internal and external cycles with external cycles discharging thermo-static liquid outside the bath raising the temperature of the external vessel.

Features

- ▶ Volume - 7 L
- ▶ Low temperature cooling pump
- ▶ Microcomputer controlled temperature
- ▶ Temperature and time LCD display
- ▶ Built-in liquid level monitor
- ▶ Refrigeration unit with relay, protector and capacitor
- ▶ Anti-corrosive circulatory system

Specifications

Model No	LRTB-A40
Volume	7 L
Temperature range	-40 °C to 200 °C
Refrigerating capacity	2329 W
Pump flow	20 L / min
Temperature stability	±0.1 °C
Head of delivery	10 m
Temperature control tier number	32
Heating power	1200