

Stainless Steel Water Distiller

LSWD-A1 series



5 L DISTILLER



10 L DISTILLER



20 L DISTILLER

Stainless Steel Water Distiller LSWD-A1 series

Stainless Steel Water distiller LSWD-A1 series are electrical floor mounted single reservoir distillation units. Stainless steel construction aids to produce germ-free distillate. Narrow piping connects heater column to heat exchanger for optimal distillate quality. Exhaust vent provides an outlet for odour or gases to be released without causing contamination.

Features

- USP grade pyrogen-free distillate
- Shunt valve regulates inlet and outlet flow
- Gas/stench exhaustion vent
- Corrosion-proof reservoir

Application

Used in laboratory, organic chemistry lab, fermentation and medical industry to produce industrial grade highly treated and distilled water suitable for laboratory usage.

Specification

MODEL NO	LSWD- A10	LSWD- A11	LSWD- A12
Capacity	5 L/hr	10 L/hr	20 L/hr
Water consumption	1:10	1:10	1:10
Configuration	Standard control	Standard control	Standard control
Water type	Regular water	Regular water	Regular water
Power	220 V ; 2.5 kW	380 V ; 2.5 kW	380 V ; 4.5 kW
	Single phase	Triplephase	Triplephase
Dimensions	365 × 368 × 845 mm	365 × 368 × 845 mm	365 × 368 × 945 mm
Weight	9 kgs	9 kgs	10kgs