

Manual Single Door
TOXIC CABINET
LTXC-B10



45 L Toxic Cabinet LTXC-B10

45 L Toxic Cabinet LTXC-B10 is a dual walled manual combustion cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 45 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LTXC-B10
Capacity	45 L
Door type	Single door manual
Adjustable shelf	1
Shelves loading capacity	50 kgs
Dimension (H x W x D)	890 x 590 x 460 mm
Weight	59 kgs