

Manual Single Door
Weak Acid and
Alkali Cabinet
Storage Cabinet
LWAC-A1 Series



Manual Single Door Weak Acid and Alkali Chemical Storage Cabinet

15 L Weak Acid and Alkali Cabinet LWAC-B10

15 L Weak Acid and Alkali Cabinet LWAC-B10 is a single walled manual Weak Acid and Alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 15 L
- Single walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-B10
Capacity	15 L
Door type	Single door manual
Adjustable shelf	1
Dimension (W x D x H)	560 x 430 x 430 mm
Weight	38 kgs

Manual Single Door Weak Acid and Alkali Chemical Storage Cabinet

45 L Weak Acid and Alkali Cabinet LWAC-B11

45 L Weak Acid and Alkali Cabinet LWAC-B11 is a single walled manual Weak Acid and Alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 45 L
- Single walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-B11
Capacity	45 L
Door type	Single door manual
Adjustable shelf	1
Dimension (W x D x H)	700 x 570 x 1055 mm
Weight	66 kgs

Manual Single Door Weak Acid and Alkali Chemical Storage Cabinet

83 L Weak Acid and Alkali Cabinet LWAC-B12

83 L Weak Acid and Alkali Cabinet LWAC-B12 is a single walled manual Weak Acid and Alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 83 L
- Single walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Three adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-B12
Capacity	83 L
Door type	Single door manual
Adjustable shelf	3
Dimension (W x D x H)	700 x 570 x 1815 mm
Weight	116 kgs