

Manual Single Door
Combustible Cabinets
LCBC-B Series



15 L Combustible Cabinet LCBC-B10

15 L Combustible Cabinet LCBC-B10 is a dual walled manual combustion cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 15 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

45 L Combustible Cabinet LCBC-B11

45 L Combustible Cabinet LCBC-B11 is a dual walled manual combustion cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 45 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specification

Model No	LCBC-B10	LCBC-B11
Capacity	15 L	45 L
Door type	Single door manual	Single door manual
Adjustable shelves	1	1
Shelves loading capacity	50 kgs	50 kgs
Dimension (H x W x D)	560 x 430 x 430 mm	890 x 590 x 460 mm
Weight	30.5 kgs	59 kgs