

# UV Pass Box

(without HEPA filters)

LUP-B1 Series



## UV Pass Box (without HEPA filters) LUP-B1 Series

---

UV Pass Box (without HEPA filters) LUP-B1 series comes with static type UV pass box equipped UV germicidal lamp, both side electrical interlocked doors, and high quality see through glass window. Pass box is made up of durable high quality 304 grade stainless steel, and updated quality parts to meet the demands of controlled clean room environment operations.

---

### Features

- Indicator lights to know the door lock/unlock condition
  - Electronically interlocked door type
  - Three on-off press buttons on the left side of box
  - Certified with CE and ISO standards
- 

### Applications

Used to transfer articles, materials between different levels of clean room across pharmaceutical, food, life science and other research, quality inspection industries.

### Specification

Model no.	LUP-B10	LUP-B11	LUP-B12
External Dimension (W*D*H)	660 x 570 x 580 mm	760 x 670 x 680 mm	860 x 770 x 780 mm
Internal Dimension (W*D*H)	500 x 500 x 500 mm	600 x 600 x 600 mm	700 x 700 x 700 mm
UV Lamp	30W*1		
Door Interlock Type	Electronical type		
Material	304 Stainless Steel		
Gross Weight	35 kg	40 kg	50 kg
Packaging Dimension (W*D*H)	700 x 700 x 780 mm	800 x 800 x 880 mm	900 x 900 x 980 mm
Power Supply	220 V, 50 Hz		